

An evaluation of monitoring possibilities of argatroban using rotational thromboelastometry and aPTT

Schött, Ulf; Rundgren, Malin; Engström, M

Published in: Journal of Thrombosis and Haemostasis

10.1111/j.1538-7836.2009.03473_2.x

2009

Link to publication

Citation for published version (APA):

Schött, U., Rundgren, M., & Engström, M. (2009). An evaluation of monitoring possibilities of argatroban using rotational thromboelastometry and aPTT. Journal of Thrombosis and Haemostasis, 7(s2), 871-872. Article PP-WE-721. https://doi.org/10.1111/j.1538-7836.2009.03473_2.x

Total number of authors:

General rights

Unless other specific re-use rights are stated the following general rights apply: Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.

 • You may not further distribute the material or use it for any profit-making activity or commercial gain

 • You may freely distribute the URL identifying the publication in the public portal

Read more about Creative commons licenses: https://creativecommons.org/licenses/

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

LUND UNIVERSITY

AN EVALUATION OF MONITORING POSSIBILITIES OF ARGATROBAN USING ROTATIONAL THROMBOELASTOMETRY AND APTT

U. Schott*1,2, M. Rundgren1, M. Engström2

¹Anesthesia and Intensive Care, University Hospital in Lund, Lund Sweden, ²Anesthesia and Intensive Care, Halmstad Central Hospital, Halmstad, Sweden


Conclusion


A significant and strong correlation between argatroban concentrations and several ROTATIONAL THROMBOELASTOMETRY (ROTEM®) parameters were found, especially so in the clinically relevant therapeutic range up to 100 s aPTT for patients with heparin induced thrombocytopenia (HIT). Rotational thrombelastometry/ thrombelastography has a potential role to increase the safety of argatroban anticoagulation in critically ill patients, like patients with heparin induced thrombocytopenia (HIT).


Introduction: The aim of the study was to studyin vitro argatroban anticoagulation with ROTEM® and activated partial thromboplastin time (aPTT).

Method: Argatroban was added in vitro to a series of citrated whole blood samples from 10 healthy volunteers to obtain whole blood concentrations of 0, 0.125, 0.25, 0.5, 1.0, 2.0, 4.0 and 8.0 mg/l. Whole blood aPTT was analysed with a Free Oscillation Rheometer (FOR or ReoRox®, Medirox Sweden) and thromboelastometry ROTEM® (Pentapharm GmbH, Munich Germany).

Results: There was a significant and strong correlation between argatroban concentrations and Clotting Time (CT in ROTEM® analysis with INTEM) (p<0.0001 and r=0.98). When we studied the correlation between aPTT and CT, we found a highly significant and strong correlation between these two analyses (p<0.0001 and r=0.97). See figs1-3.







References

Shore-Lesserson, Evidence based coagulation monitors: heparin monitoring, thromboelastography, and platelet function, Semin Cardiothorac Vasc Anesth 9 (2005), pp. 41–52.

Taketomi T, Szlam F, Vinten-Johansen, J, Levy JH, Tanaka KA. Thrombin-activated thrombelastography for evaluation of thrombin interaction with thrombin inhibitors. Blood Coagulation & Fibrinolysis: Volume 18(8) December 2007p 761-76.

Sucker C, ZotZ B, Görlinger K, Hartmann M. Rotational thrombelastography for bedside monitoring of recombinant hirudin. Acta Anaesthesiol Scand 2008;52: 358-362.

Kouerinis IA, El-Ali M, Theakos N, Dedeillas P. Can thrombelstography predict which patients with heparin-induced thrombocytopenia may suffer thrombotic complications of type II. Eur J Cardio-Thoracic Surg 2007;32: 544-546.